

PRODUCT INFORMATION

- **Complete CPAP Oxygen Delivery System**
- **Simple to use**
- **Reduces Oxygen Tank Consumption**
- **Reduces the Need for Intubation in CHF Patients**
- **Reduces the Cost of Care**



PORTOVENT
Product Features

Superior Technology:

- Non-invasive Positive Pressure Ventilation
- Totally Portable
- No PEEP valve
- Completely self contained design
- Adjustable Pressure Settings with Gauge
- Internal Manometer Displays Real Time Airway Pressure

Accuracy:

- Automatic leak compensation delivers constant CPAP
- Comfort Seal mask provides superior seal and comfort for the patient

Always Ready:

- No recharging required
- No batteries to maintain
- No electric power required in an emergency
- Totally Pneumatic

Regulatory Clearances

- U.S. FDA 510(k)
- CE-Mark
- UL/CSA Classification

Soft Seal Masks

Advanced Design (Patent Pending):

- Uniquely formulated, latex free, non-reticulating foam seal
 - Virtually leak free and conforms to a variety of facial morphologies, beards and tubing.
 - Very comfortable for the patient

Patient Exhalation Valve (Patent Pending):

- Reduces the inspiratory work of breathing while maintaining CPAP during exhalation
- Large 360° "swivel" port directs exhalation and can accommodate a standard 30mm bacterial/viral filter

PORTOVENT

CPAP Oxygen Delivery System

EMT APPLICATIONS

Lowest Cost, Reduces Risks, Improves Outcomes

The Portable EMT CPAP Oxygen Delivery System

- ✓ Reduces Amount of Oxygen Consumed
- ✓ Less Traumatic, Non-Invasive CPAP Technology
- ✓ Reduces Hospital Length of Stays and Risks of Infection for Many Patients

EMERGENCY DEPARTMENT APPLICATIONS

Lowest Cost, Reduces Risks, Improves Outcomes

The Portable ED Oxygen Delivery System

- ✓ Non-Invasive CPAP Technology Reduces The Need For High Cost Ventilators Who Only Require Positive Pressure Ventilation
- ✓ Clinically Proven to Lower the Work of Breathing
- ✓ Reduces Hospital Length of Stays and Risk of Infection for Many Patients
- ✓ Improves Quality of Patient Care



Note: All features are nominal and subject to change without notice.

PortO₂Vent is a trademark of Emergent Respiratory Products.

PortO₂Vent Rev. 1.0
©2004 Emergent. All Rights Reserved

**For More Information About
Emergent Respiratory Products
Visit our website @
WWW.ERESP.COM**